

Name _____ Date _____ Partners _____

EQUIPOTENTIAL AND THE ELECTRIC FIELD

Statement of Purpose

The purpose of the following experiment is the observation of an electric field by quantitative measurement. By determining the voltage at various points between two electrodes, we are able to map equipotential lines; these lines can be used to map the electric field lines, which are at all points perpendicular to the equipotential.

Notes on the Multimeter

Source: Philips, W.D. "2P: Using a Multimeter." *Design Electronics*. DOCTRONICS Educational Publications. 1999. <http://www.doctrronics.co.uk/meter.htm>.

The multimeter we will be using is the Wavetek Meterman 27XT Digital Multimeter. The multimeter is an extremely useful tool when analyzing electrical circuits because it is able to operate in three distinct modes: by measuring current (as an ammeter), by measuring voltage (as a voltmeter), and by measuring resistance (as an ohmmeter). When the multimeter is functioning as an ammeter, the current passing through it must be identical to the current in the circuit loop for which the measurement is being taken; therefore, the ammeter/multimeter must be connected in series. This function requires a low resistance to avoid modifying the circuit. When the multimeter is functioning as a voltmeter, the voltage drop across the multimeter must be equivalent to that across the circuit branch being examined; therefore, the voltmeter/multimeter must be connected in parallel, as the voltage drop across all branches in a circuit is uniform. This function requires a high resistance to avoid shorting the circuit. The multimeter may additionally function as an ohmmeter when connected across a circuit element which is separate from the circuit. For the purpose of this experiment, we will use the multimeter as a voltmeter.

Experimental Setup

Equipment	Qty
HP 6216A DC Power Supply	1
Wavetek Meterman 27XT Multimeter	1
Brass cylindrical "point" electrodes	2
Aluminum plate electrodes	2
Base and support rod	2
Three finger clamp	2
Stacking banana plugs (black)	2
Banana plugs (one red, one black)	2
Square plastic dish, 8x8 (in.)	1
Waterproofed grid sheet	1
Duplicate grid sheet (not waterproof)	3

Colored markers	5
Water	

I. Two Point Electrodes

1. With the waterproofed grid sheet placed in the bottom of the plastic dish, fill the dish with water to a depth of one to two centimeters. You may want to tape the grid down so it will not slide around as you take data.
2. Turn the Wavetek Meterman Multimeter knob to 20V and select DC mode.
3. With the three-fingered prongs attached to the base and support rods, place the two “point” electrodes in the clamps. Position the holders so that the electrodes are roughly equidistant from their respective sides of the dish. The lower end of each electrode should be submerged in the water. Use a red banana plug to connect one electrode to the positive terminal of the DC Power Supply and a stacking banana plug to connect the other electrode to the negative terminal.
4. Using a second stacking banana plug, connect the black electrode to the COM input of the Wavetek Meterman Multimeter.
5. Plug the multimeter probe lead into the V Ω input of the Wavetek Meterman Multimeter.




Figure 1

6. **Potential Difference:** Turn on the DC source. Set the input to 5 Volts. This should supply a current and create a potential difference across the electrodes. Bring the probe in contact under water with one of the two electrodes and check the voltage reading with the multimeter.

Potential (1st Electrode): _____

Question 1: What voltage reading do you expect to find for the other electrode?

Use the probe to check your prediction.

Potential (2nd Electrode): _____

7. **Central Equipotential Line:** Locate a point midway between the two electrodes (as close as possible to the exact center) and measure the potential here by touching the probe to the bottom of the container. Keep in mind that the voltage reading may not stabilize immediately, so you may need to hold the probe still for several seconds before getting an accurate reading. Choose a color to represent this voltage and record the voltage and corresponding color in the table below (Table 1); also record the grid position. Use the marker to mark the position of this first potential on the duplicate grid, which you should label "Grid 1." Now use the probe and multimeter to locate a second point with the same potential. Write this grid position in Table 1 and use the marker to record it on the duplicate grid. Repeat this process until you have located several points on each side of the central point (i.e., the point on the axis of the electrodes). List each grid position in the Table and on the duplicate grid. A line connected through these points should represent your central equipotential line; draw this line on the duplicate grid.

Color	Voltage (V)	Grid Positions

Table 1

8. **Off-Center Equipotential Lines:** Repeat the process described in Step 7, but this time chose an initial voltage point between the electrodes on either side of the center point. Again, locate several points on either side of the electrode axis and record these values in different color ink than that used in Step 7. Remember to record the color, voltage, and grid positions in Table 1.

Prediction 1: What would a connecting line between these points resemble? Sketch your prediction below.

Repeat this process for all five colors, twice on each side of the central equipotential line.

9. Field Lines: Electric field lines are perpendicular to equipotential lines at all points.

Prediction 2: Before drawing the field lines on your data sheet, predict the field line configuration by drawing it below:

Now draw short dashes on the duplicate grid perpendicular to your equipotential lines at several points along each line. Connect these dashes in order to draw electric field lines.

II. Point and Plane Electrodes

1. With the power supply turned off to avoid electrical shock, replace one of the “point” electrodes with a plane electrode. The plane electrode balances on its narrow side on the bottom of the dish, so the base and rod are not needed. Insert the banana plug formerly connected to the “point” electrode into the hole in the corner of the plane electrode.
2. Turn the voltage back on and repeat the procedure outlined above for this configuration. The central electrode axis should now be located between the electrode and the center of the plate.

Prediction 3: How do you expect the equipotential and the field lines to appear with this electrode configuration? Sketch your predictions below.

Test your predictions by following the process described above. Use a new duplicate grid (label it “Grid 2”) and record data in Table 2.

Color	Voltage (V)	Grid Positions

Table 2

III. Two Plane Electrodes

1. With the power supply turned off to avoid electrical shock, replace the remaining “point” electrode with a second plane electrode. The apparatus should now contain two parallel planes within the configuration described in the initial laboratory setup.
2. Turn the voltage on and repeat the procedure outlined in Part I for this configuration.

Prediction 4: How do you expect the equipotential and the field lines to appear with these electrode shapes? Sketch your predictions below.

Test your predictions. Use a new duplicate grid (label it "Grid 3") and record data in Table 3.

Color	Voltage (V)	Grid Positions

Table 3

Grid sheet

	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46
47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69
70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92
93	94	95	96	97	98	99	100	101	102	103	104	105	106	107	108	109	110	111	112	113	114	115
116	117	118	119	120	121	122	123	124	125	126	127	128	129	130	131	132	133	134	135	136	137	138
139	140	141	142	143	144	145	146	147	148	149	150	151	152	153	154	155	156	157	158	159	160	161
162	163	164	165	166	167	168	169	170	171	172	173	174	175	176	177	178	179	180	181	182	183	184
185	186	187	188	189	190	191	192	193	194	195	196	197	198	199	200	201	202	203	204	205	206	207
208	209	210	211	212	213	214	215	216	217	218	219	220	221	222	223	224	225	226	227	228	229	230
231	232	233	234	235	236	237	238	239	240	241	242	243	244	245	246	247	248	249	250	251	252	253
254	255	256	257	258	259	260	261	262	263	264	265	266	267	268	269	270	271	272	273	274	275	276
277	278	279	280	281	282	283	284	285	286	287	288	289	290	291	292	293	294	295	296	297	298	299
300	301	302	303	304	305	306	307	308	309	310	311	312	313	314	315	316	317	318	319	320	321	322
323	324	325	326	327	328	329	330	331	332	333	334	335	336	337	338	339	340	341	342	343	344	345
346	347	348	349	350	351	352	353	354	355	356	357	358	359	360	361	362	363	364	365	366	367	368
369	370	371	372	373	374	375	376	377	378	379	380	381	382	383	384	385	386	387	388	389	390	391
392	393	394	395	396	397	398	399	400	401	402	403	404	405	406	407	408	409	410	411	412	413	414
415	416	417	418	419	420	421	422	423	424	425	426	427	428	429	430	431	432	433	434	435	436	437
438	439	440	441	442	443	444	445	446	447	448	449	450	451	452	453	454	455	456	457	458	459	460
461	462	463	464	465	466	467	468	469	470	471	472	473	474	475	476	477	478	479	480	481	482	483
484	485	486	487	488	489	490	491	492	493	494	495	496	497	498	499	500	501	502	503	504	505	506
507	508	509	510	511	512	513	514	515	516	517	518	519	520	521	522	523	524	525	526	527	528	529
530	531	532	533	534	535	536	537	538	539	540	541	542	543	544	545	546	547	548	549	550	551	552
553	554	555	556	557	558	559	560	561	562	563	564	565	566	567	568	569	570	571	572	573	574	575
576	577	578	579	580	581	582	583	584	585	586	587	588	589	590	591	592	593	594	595	596	597	598
599	600	601	602	603	604	605	606	607	608	609	610	611	612	613	614	615	616	617	618	619	620	621
622	623	624	625	626	627	628	629	630	631	632	633	634	635	636	637	638	639	640	641	642	643	644
645	646	647	648	649	650	651	652	653	654	655	656	657	658	659	660	661	662	663	664	665	666	667
668	669	670	671	672	673	674	675	676	677	678	679	680	681	682	683	684	685	686	687	688	689	690
691	692	693	694	695	696	697	698	699	700	701	702	703	704	705	706	707	708	709	710	711	712	713
714	715	716	717	718	719	720	721	722	723	724	725	726	727	728	729	730	731	732	733	734	735	736
737	738	739	740	741	742	743	744	745	746	747	748	749	750	751	752	753	754	755	756	757	758	759
760	761	762	763	764	765	766	767	768	769	770	771	772	773	774	775	776	777	778	779	780	781	782
783	784	785	786	787	788	789	790	791	792	793	794	795	796	797	798	799	800	801	802	803	804	805
806	807	808	809	810	811	812	813	814	815	816	817	818	819	820	821	822	823	824	825	826	827	828
829	830	831	832	833	834	835	836	837	838	839	840	841	842	843	844	845	846	847	848	849	850	851
852	853	854	855	856	857	858	859	860	861	862	863	864	865	866	867	868	869	870	871	872	873	874
875	876	877	878	879	880	881	882	883	884	885	886	887	888	889	890	891	892	893	894	895	896	897
898	899	900	901	902	903	904	905	906	907	908	909	910	911	912	913	914	915	916	917	918	919	920
921	922	923	924	925	926	927	928	929	930	931	932	933	934	935	936	937	938	939	940	941	942	943
944	945	946	947	948	949	950	951	952	953	954	955	956	957	958	959	960	961	962	963	964	965	966
967	968	969	970	971	972	973	974	975	976	977	978	979	980	981	982	983	984	985	986	987	988	989
990	991	992	993	994	995	996	997	998	999	1000	1001	1002	1003	1004	1005	1006	1007	1008	1009	1010	1011	1012
1013	1014	1015	1016	1017	1018	1019	1020	1021	1022	1023	1024	1025	1026	1027	1028	1029	1030	1031	1032	1033	1034	1035
1036	1037	1038	1039	1040	1041	1042	1043	1044	1045	1046	1047	1048	1049	1050	1051	1052	1053	1054	1055	1056	1057	1058
1059	1060	1061	1062	1063	1064	1065	1066	1067	1068	1069	1070	1071	1072	1073	1074	1075	1076	1077	1078	1079	1080	1081
1082	1083	1084	1085	1086	1087	1088	1089	1090	1091	1092	1093	1094	1095	1096	1097	1098	1099	1100	1101	1102	1103	1104
1105	1106	1107	1108	1109	1110	1111	1112	1113	1114	1115	1116	1117	1118	1119	1120	1121	1122	1123	1124	1125	1126	1127
1128	1129	1130	1131	1132	1133	1134	1135	1136	1137	1138	1139	1140	1141	1142	1143	1144	1145	1146	1147	1148	1149	1150
1151	1152	1153	1154	1155	1156	1157	1158	1159	1160	1161	1162	1163	1164	1165	1166	1167	1168	1169	1170	1171	1172	1173
1174	1175	1176	1177	1178	1179	1180	1181	1182	1183	1184	1185	1186	1187	1188	1189	1190	1191	1192	1193	1194	1195	1196
1197	1198	1199	1200	1201	1202	1203	1204	1205	1206	1207	1208	1209	1210	1211	1212	1213	1214	1215	1216	1217	1218	1219

